**PURPOSE**

This unit enables to inhale a breathing mixture which is different from the one constituting the ambient atmosphere.

It has been specially designed to breathe oxygen-rich mixtures used in hyperbaric oxygen therapy or by divers during chamber decompression.

Masks available in 3 sizes : S, M and L.

DESCRIPTION

The breathing set consists of an oral-nasal mask equipped with a Demand Regulator and an exhaust valve.

The regulator balanced demand-valve can supply a large gas-flow for a minimum inhaling effort.

In order to prevent oxygen build-up in the atmosphere of the chamber, thus minimising fire risks, the Demand Regulator Breathing Mask should be connected to an Oxygen Dump Valve, such as the LB78 Oxygen Dump Valve ref. 421 00211.

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CHARACTERISTICS

- Individual breathing equipment
- Fully rated for use with pure oxygen under pressure
- Gas Supply relative pressure : 3 to 8 bar
- Valve opening (overpressure) : 3 millibar
- Maximum instant flow rate : > 400 l/min
for a 5 millibar depression
- Demand regulator Inlet : 9/16" SAE
M12x100 adaptor available
- Mask outlet is a dia. 22 mm straight hose connector
that may easily be adapted to any LP hoses
- Housing made of medical POM-C and Stainless Steel
- Mechanism of chrome plated brass and Stainless Steel
- Oro-nasal Mask and Diaphragm made of silicone
- O' ring sealing
- Simple and rugged, no adjustment needed
- Weight : 450 g

Corrugated Hose (sold separately):

- Heavy-duty LP Hose to connect the Demand Regulator
Breathing Mask to the Dump valve
- Made of Medical Silicone
- Length : 800 mm

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
MAINTENANCE

- Thanks to its design, the demand regulator requires no adjustment.
- For optimal functioning, it should be cleaned every 6 months and the O'rings should be slightly lubricated with an inflammable liquid grease.
- O'rings, diaphragms (regulator and mask), the vulcanised seat and the nozzle should be changed every year (maintenance kit 411 01027).

DISMANTLING

- Unscrew 260 09734 to remove the clamp 411 09643.
- Take off the cover 411 09636.
- Take off the cone 411 09639.
- Take out diaphragm 411 09640.
- Unscrew the seat holder 411 09635.
- Unscrew the vulcanised seat 411 09691 from the seat holder 411 09635.
- Slowly unscrew the plug 411 09739 to accommodate extension of the spring 400 11886.
- Remove the clip 411 09633 in order to remove the nozzle 411 09634 from the connection pipe 411 09631.
- Disengage lever 411 09632 sliding it in the groove of the connection pipe 411 09631.
- Remove the connection pipe 411 09631.

CLEANING

- Clean all parts in soapy water.
- Rinse thoroughly. Dry using clean and fresh air.
- Lubricate slightly nozzle 411 09634 and O'rings using an inflammable liquid grease.

ASSEMBLING AND SETTING

- Assembling in reverse order to the dismantling.
- At the end of the assembling, screw the seat holder 411 09635 in the connection pipe 411 09631 until the lever 411 09632 detaches from the demand regulator casing 411 11775 (get a movement of +/- 1 mm at the tip of the lever).

In the event of a leak,

- Check movement of +/- 1 mm at the tip of the lever.
- Check all O'rings.
- Check vulcanised seat 411 09691.
- Check stainless steel nozzle 411 09634, inspect with particular care the edge which shall not be scratched.
- Check surface condition of the connection pipe 411 09631 where the nozzle 411 09634 translates.

The demand regulator hard to blow through,

- Clean and lubricate it.
- Check movement of +/- 1 mm at the tip of the lever.
- Check diaphragm 411 09640 which may have lost its flexibility, in which case, it should be changed.
- Check spring 400 11886.

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